

Gentec[®] offers a complete line of Oxygen and Air Medical Gas Flowmeters to meet your clinical needs



FM187B-15L-00



FM197B-15L-DH-P

Oxygen Flowmeter w/ Power Take-off
& Chemetron Adapter

- Gentec flowmeters have a 0-15 LPM flow range, with a minimum of 65 LPM flush flow rate
- The bases are available in either nickel plated solid brass or aluminum
- Durable polycarbonate hood and metering tube
- Stainless steel ball float
- Flowmeters are all calibrated at 50 psig and 70° F
- Flowmeters are all back-pressured compensated
- Color-coded for air and oxygen
- MRI (optional)
- Dual flowmeter (available)
- ISO flowmeters air/oxygen (available)
- 5-year warranty (see warranty statement for details)

Please refer to the ordering information on the back side of the page, or contact your GENTEC[®] distributor for assistance with selecting your new GENTEC[®] flowmeters

Flowmeter Ordering Information

| FLOWMETER TYPE ADAPTER | AIR FLOWMETER 0-15 LPM | AIR FLOWMETER W/ POWER TAKE-OFF | AIR DUAL FLOWMETER 0-15 LPM | OXYGEN FLOWMETER 0-15 LPM | OXYGEN W/ POWER TAKE-OFF | OXYGEN DUAL FLOWMETER 0-15 LPM |
|-------------------------------|-----------------------------------|------------------------------------|-----------------------------------|---------------------------------|--------------------------------|--------------------------------------|
| | STANDARD INLET 1/8" NPT FEMALE | | | | | |
| | FM187B-15L-00 | FM187B-15L-00-P | FMD187B-15L-00 | FM197B-15L-00 | FM197B-15L-00-P | FMD197B-15L-00 |
| | FM187B-15L-OH | FM187B-15L-OH-P | FMD187B-15L-OH | FM197B-15L-OH | FM197B-15L-OH-P | FMD197B-15L-OH |
| | FM187B-15L-CH | FM187B-15L-CH-P | FMD187B-15L-CH | FM197B-15L-CH | FM197B-15L-CH-P | FMD197B-15L-CH |
| | FM187B-15L-DH | FM187B-15L-DH-P | FMD187B-15L-DH | FM197B-15L-DH | FM197B-15L-DH-P | FMD197B-15L-DH |
| | FM187B-15L-DS | FM187B-15L-DS-P | FMD187B-15L-DS | FM197B-15L-DS | FM197B-15L-DS-P | FMD197B-15L-DS |
| | FM187B-15L-GS | FM187B-15L-GS-P | FMD187B-15L-GS | FM197B-15L-GS | FM197B-15L-GS-P | FMD197B-15L-GS |
| | FM187B-15L-BS | FM187B-15L-BS-P | FMD187B-15L-BS | FM197B-15L-BS | FM197B-15L-BS-P | FMD197B-15L-BS |
| | FM187B-15L-FS | FM187B-15L-FS-P | FMD187B-15L-FS | FM197B-15L-FS | FM197B-15L-FS-P | FMD197B-15L-FS |

Options:

- For aluminum flowmeters, change the B to an A (example: FM197B-15L-OH to FM197A-15L-OH)
- For dual flowmeters with 2 power take-offs, add -2P (example: FMD197B-15L-OH-2P)
- For ISO flowmeters add an E (example: FM197A-15L-OHE)
- For MRI flowmeters add MRI to the end (example: FM197A-15L-OH-MRI)

Specifications

| | |
|----------------------------|---|
| Weight (no adapter) | Aluminum Body: 4.24 oz/122 gm Brass Body: 6.64 oz/188 gm |
| Increment: | 0-6 LPM in ½ LPM increment, 6-15 LPM in 1 LPM Increment |
| Flush Flow Rate: | 65 LPM (depending on supply flow) |
| Max Inlet Pressure: | 100 psig |
| PTO Flow Rate: | Over 100 LPM (depending on supply flow) |