


## Features

- The Ceiling Retractable Power Column is designed and manufactured according to IEC 60601-1 Medical Electrical Equipment.
- Our equipment are UL Listed for U.S. and Canada.
- All parts are fabricated using CAD/CAM software and CNC machinery to ensure the precision and quality of all the assembly. Our company is ISO 9001:2015.
- Each part not fabricated by GENTEC® are UL Listed.



## Body Specifications

- Exterior: Stainless steel 304 GA18
- Internal structure: Galvanized Steel CA18
- Finish: Satin
- Weight: 150 pounds
- Internal division between gas and electrical
- Ceiling mounting
- Service access: Back and bottom



## Technical Specifications

- Gas pressure: Oxygen: 60 PSI, Med Air: 60 PSI, Med Vac: inHg
- Type of gas outlet: DISS
- Connection to the gas pipeline: Conductive hose 1/4" NFWA color code, DISS adaptor.
- Electrical: 115/220 VAC - 3 circuits
- Internal wiring gauge: 12 AWG
- Type of wire: HF - Halogen free, FR - Flame Retardant, LS - Low Smoke
- Electrical outlets: Hospital grade, 20 A, NEMA 5-20-R, 2 pole, 3 wires

## Structure Required

