


ART WALL ENCLOSURE

Features

- Designed for hospital rooms where design is emphasized, replaces traditional bedhead units
- The painting sits on a sliding mechanism which can hide the medical equipment while they are not in use
- The art piece and exterior design can be customized to meet or match a variety of design requirements
- A wide range of medical products & accessories can be custom-ordered
- Gas and electric channels are insulated to ensure safety



Dimensions



Note: All dimensions are reference.

Ordering Information

3721 - P01 - GS		
Series	Canvas Style	Medical Gas Outlet
3721	P01: Portrait	OH: Ohmeda® adapter
	P02: Landscape	DS: DISS hex adapter
	P03: Abstract	DH: DISS hand tight adapter
		CH: Chemetron® adapter
		PB: Puritan-bennett® adapter
		FS: French (NF S 90-116) adapter
		GS: German (DIN 13260-2) adapter
		BS: British (BS5682-1998) adapter
		JIS: Japanese Style adapter
		SIS: Australian (AS2896) adapter

Please contact Genstar for more ordering information.