



Product Highlights

- Portable, Lightweight & Rugged
- Adjustable suction levels (300 mmHg or 600 mHg)
- One-hand operation
- No electric power required
- 300 ml disposable collection canister
- Simple to use
- Adult & infant suction catheters included
- 5 year warranty

Kit Includes:

- 1 disposable canister
- 2 suction catheters (6mm and 12mm)
- 1 reusable pump

Model Number	Description
MSP-1	Handi Pump Kit= includes reusable pump, one canister, adult and child catheter
MSP-1-PK	Replacement canister with 2 suction catheters (6mm and 12mm)

Product Specifications			
Maximum vacuum	600 mmHg	Suction Levels	300mmHg & 600mmHg
Flow rate (maximum)	20L/min	Operating Range	-4° F to 122° F
Collection canister connection	17 mm	Storage Range	-40° F to 140° F
Collection canister volume	300 ml	Size	185 x 64 x 168mm
Catheter connection OD/ID	13 mm/6 mm	Weight	230g (8.11 oz)