



























Medical Hoses & Adapters Ordering Information (International Edition)



CATA#: MED17V1509
© Gentec (Shanghai) Corporation. September 2015.

Medical Gas Adapters

Ohmeda® Adapters














| Ohmeda® |  |  |  |  |  |
|-------------------------|--|--|--|--|--|
| USA Standard | Oxygen | Vacuum | Medical Air | Nitrous Oxide | WAGD |
| 1/8" NPT |  204301610 3111U-OXY-NT2 |  204301620 3111U-VAC-NT2 |  204301630 3111U-AIR-NT2 |  204301640 3111E-N2O-NT2 |  204301650 3111U-WAG-NT2 |
| 1/4" or 5/16" Hose Barb |  204301710 3111U-OXY-HB4 |  204301725 3111U-VAC-HB5 |  204301730 3111U-AIR-HB4 |  204301740 3111U-N2O-HB4 |  204301750 3111U-WAG-HB4 |
| ISO Standard | Oxygen | Vacuum | Medical Air | Nitrous Oxide | WAGD |
| 1/8" NPT |  204301611 3111E-OXY-NT2 |  204301621 3111E-VAC-NT2 |  204301631 3111E-AIR-NT2 |  204301640 3111E-N2O-NT2 | |
| 1/4" or 5/16" Hose Barb |  204301711 3111E-OXY-HB4 |  204301726 3111E-VAC-HB5 |  204301731 3111E-AIR-HB4 |  204301740 3111U-N2O-HB4 |  204301756 3111E-WAG-HB5 |

Chemetron® Adapters

























| Chemetron® |  |  |  |  |  |
|--------------|---|---|---|--|--|
| USA Standard | Oxygen | Vacuum | Medical Air | Nitrous Oxide | WAGD |
| 1/8" NPT |  201160040 3103-OXY-NT2 |  201160050 3103-VAC-NT2 |  201160020 3103-AIR-NT2 |  201160030 3103-N2O-NT2 |  201160110 3103-WAG-NT2 |
| ISO Standard | Oxygen | Vacuum | Medical Air | | |
| 1/8" NPT |  204303011 3103E-OXY-NT2 |  204303021 3103E-VAC-NT2 |  204303031 3103E-AIR-NT2 |  201160030 3103-N2O-NT2 | |

Note: If the panel is too large, you cannot directly use the machine lock. The HB need use NT2 + FNT2-HB.

Puritan-Bennett® Adapters

























| Puritan Bennett® |  |  |  |  |  |
|-------------------------|---|---|---|---|---|
| | Oxygen | Vacuum | Medical Air | Nitrous Oxide | WAGD |
| 1/8" NPT |  |  |  |  |  |
| Chrome-plated Brass | 301160290 3104-OXY-NT2 | 301160286 3104-VAC-NT2 | 301160260 3104-AIR-NT2 | 301160250 3104-N2O-NT2 | 301160270 3104-WAG-NT2 |
| Stainless Steel | 301160910 3114-OXY-NT2 | 301160920 3114-VAC-NT2 | 204311430 3114-AIR-NT2 | 301161351 3114-N2O-NT2 | / |
| 1/4" or 5/16" Hose Barb |  |  |  |  |  |
| Chrome-plated Brass | 301160580 3104-OXY-HB4 | 301160560 3104-VAC-HB5 | 301160530 3104-AIR-HB4 | 301160570 3104-N2O-HB4 | 301160550 3104-WAG-HB5 |
| Stainless Steel | C53300410 3114-OXY-HB4 | C53300420 3114-VAC-HB5 | C53300431 3114-AIR-HB4 | C53300440 3114-N2O-HB4 | C53300451 3114-WAGD-HB5 |

DISS Hex Adapters

























| DISS Hex |  |  |  |  |
|-------------------------|---|---|---|---|
| | Oxygen | Vacuum | Medical Air | Nitrous Oxide |
| 1/8" NPT |  |  |  |  |
| | 204302110 3102-OXY-NT2 | 204302120 3102-VAC-NT2 | 204302130 3102-AIR-NT2 | 204302140 3102-N2O-NT2 |
| 1/4" or 5/16" Hose Barb |  |  |  |  |
| | 204302612 3102L-OXY-HB4 | 204302622 3102L-VAC-HB5 | 204302632 3102L-AIR-HB4 | 204302642 3102L-N2O-HB4 |
| DISS Hex |  |  |  |  |
| | WAGD | Carbon Dioxide | Nitrogen | Instrument Air |
| 1/8" NPT |  |  |  |  |
| | 204302150 3102-WAG-NT2 | 204302160 3102-CO2-NT2 | 204302170 3102-NIT-NT2 | 204302180 3102-IA-NT2 |
| 1/4" or 5/16" Hose Barb |  |  |  |  |
| | 204302650 3102L-WAGD-HB5 | 204302660 3102L-CO2-HB4 | 204302672 3102L-NIT-HB4 | 204302682 3102L-IA-HB4 |

Medical Gas Adapters







DISS Manual Adapters (USA Standard)

| | | | | |
|-------------------------|--|---|--|---|
| DISS Manual |  |  |  |  |
| USA Standard | Oxygen | Vacuum | Medical Air | Nitrous Oxide |
| 1/8" NPT |  204312110 3222-OXY-NT2 |  204312122 3222U-VAC-NT2 |  204312132 3222U-AIR-NT2 |  204312141 3222E-AIR-NT2 |
| 1/4" or 5/16" Hose Barb |  204312210 3222-OXY-HB4 |  204312224 3222U-VAC-HB5(wh) |  204312232 3222U-AIR-HB4 |  204312241 3222E-N2O-HB4 |
| DISS Manual |  |  |  |  |
| USA Standard | WAGD | Carbon Dioxide | Nitrogen | Instrument Air |
| 1/8" NPT |  204312150 3222-WAG-NT2 |  204312161 3222E-CO2-NT2 |  3222E-NIT-NT2 |  204322180 3222-IA-NT2 |
| 1/4" or 5/16" Hose Barb |  204312255 3222-WAG-HB5 |  204312261 3222E-CO2-HB4 |  204312271 3222E-NIT-HB4 |  204312280 3222-IA-HB4 |











DISS Manual Adapters (ISO Standard)

| | | | | |
|-------------------------|---|---|--|---|
| DISS Manual |  |  |  |  |
| ISO Standard | Oxygen | Vacuum | Medical Air | Nitrous Oxide |
| 1/8" NPT |  204312111 3222E-OXY-NT2 |  204312121 3222E-VAC-NT2 |  204312131 3222E-AIR-NT2 |  204312141 3222E-AIR-NT2 |
| 1/4" or 5/16" Hose Barb |  204312211 3222E-OXY-HB4 |  204312226 3222E-VAC-HB5 |  204312231 3222E-AIR-HB4 |  204312241 3222E-N2O-HB4 |
| DISS Manual |  |  |  |  |
| ISO Standard | WAGD | Carbon Dioxide | Nitrogen | Instrument Air |
| 1/8" NPT |  204312151 3222E-WAG-NT2 |  204312161 3222E-CO2-NT2 |  3222E-NIT-NT2 |  204322180 3222-IA-NT2 |
| 1/4" or 5/16" Hose Barb |  204312256 3222E-WAG-HB5 |  204312261 3222E-CO2-HB4 |  204312271 3222E-NIT-HB4 |  204312280 3222-IA-HB4 |













DISS Adapter with Check Valve

| DISS | Oxygen | Vacuum | Medical Air | Nitrous Oxide |
|------------------------------|--|--|---|--|
| 1/8" NPT with Check Valve |  204031611 DISS-O-NT2-CHK |  204031620 DIDD-V-1/8MNPT |  204031633 DISS-A-NT2/chk |  204031673 DISS-N-NT2/chk |
| 1/4" or 5/16" Hose Barb |  301163511 DISS-O-HB4 |  301163521 DISS-V-HB5 |  301163531 DISS-A-HB4 |  301163541 DISS-2-HB4 |













FS Adapters

| FS |  |  |  |  |
|----------------------------|---|---|--|---|
| | Oxygen | Vacuum | Medical Air | Nitrous Oxide |
| 1/8" NPT |  204306052 3105E-O-NT2 |  204306062 3105E-V-NT2 |  204306072 3105E-A-NT2 |  204306082 3105E-2-NT2 |
| 1/4" or 5/16" Hose Barb |  204306051 3105E-O-HB4 |  204306061 3105E-V-HB5 |  204306071 3105E-A-HB4 |  204306081 3105E-2-HB4 |





















SIS Adapters

| SIS |  |  |  |  |
|----------------------------|--|--|--|--|
| | Oxygen | Vacuum | Medical Air | Nitrous Oxide |
| 1/8" NPT |  204310911 3109-OXY-NT2 |  204310921 3109-VAC-NT2 |  204310931 3109-AIR-NT2 |  204310941 3109-N2O-NT2 |
| 1/4" or 5/16" Hose Barb |  204305610 3109-OXY-HB4 |  204305620 3109-VAC-HB5 |  204305630 3109-AIR-HB4 |  204305640 3109-N2O-HB4 |

JIS Adapters






















| JIS |  |  |  |  |
|-------------------------|---|---|--|---|
| | Oxygen | Vacuum | Medical Air | Nitrous Oxide |
| 1/8" NPT |  204310810 3108-OXY-NT2 |  204310820 3108-VAC-NT2 |  204310830 3108-AIR-NT2 |  204310840 3108-N2O-NT2 |
| 1/4" or 5/16" Hose Barb |  204311811 3108-OXY-HB4H |  204311821 3108-VAC-HB5H |  204311831 3108-AIR-HB4H |  204311841 3108-N2O-HB4H |

GS Adapters

| GS |  |  |  |  |  |  |  |
|-------------------------|--|--|--|--|--|--|---|
| | Oxygen | Vacuum | Medical Air | Nitrous Oxide | Carbon Dioxide | Nitrogen | Instrument Air |
| 1/4" or 5/16" Hose Barb |  305631010 3106-OXY-HB4 |  305631020 3106-VAC-HB5 |  305631030 3106-AIR-HB4 |  305631040 3106-N2O-HB4 |  305631060 3106-CO2-HB4 |  305631070 3106-N2-HB4 | |
| 1/4" or 5/16" Hose Barb |  204311610 3116-OXY-HB4H |  204311625 3116-VAC-HB5H |  204311630 3116-AIR-HB4H |  204311640 3116-N2O-HB4H |  204311660 3116-CO2-HB4H |  204311670 3116-NIT-HB4H |  204311680 3116-AIR8-HB4H |













Note: The material of the 3106 series is chrome-plated brass.

BS Adapters
















| BS |  |  |  |  |  |  |  |
|-------------------------|--|--|--|--|---|--|--|
| | Oxygen | Vacuum | Medical Air | Nitrous Oxide | Carbon Dioxide | Nitrogen | Instrument Air |
| 1/8" NPT |  |  |  |  |  |  |  |
| Chrome-plated Brass | 301461140 3107-OXY-NT2 | 301461160 3107-VAC-NT2 | 301461170 3107-AIR-NT2 | 301461150 3107-N2O-NT2 | / | / | / |
| Stainless Steel | 204311713 3117-OXY-NT2 | 204311722 3117-VAC-NT2 | 204311731 3117-AIR-NT2 | 204311741 3117-N2O-NT2 | 204311761 3117-CO2-NT2 | 204311781 3117-N2-NT2 | 204311771 3117-MA7-NT2 |
| 1/4" or 5/16" Hose Barb |  |  |  |  |  |  |  |
| Chrome-plated Brass | 301461100 3107-OXY-HB4 | 301461120 3107-VAC-HB5 | 301461130 3107-AIR-HB4 | 301461110 3107-N2O-HB4 | / | / | 301461180 3107-MA7-HB4 |
| Stainless Steel | 204311715 3117-OXY-HB4 | 204311725 3117-VAC-HB5 | 204311735 3117-AIR-HB4 | 204311745 3117-N2O-HB4 | 204311760 3117-CO2-HB4 | 204311780 3117-N2-HB4 | 204311770 3117-MA7-HB4 |

Couplers

Ohmeda® Couplers













| Ohmeda® |  |  |  |  |
|-------------------------|---|---|--|---|
| USA Standard | Oxygen | Vacuum | Medical Air | Nitrous Oxide |
| 1/8" NPT |  A65401012 3911U-O-NT2 |  A65401022 3911U-V-NT2 |  A65401132 3911U-A-NT2 |  3911U-2-NT2 |
| 1/4" or 5/16" Hose Barb |  A65401011 3911U-O-HB4 |  A65401025 3911U-V-HB5 |  A65401131 3911U-A-HB4 |  A65401041 3911U-2-HB4 |

DIN Couplers













| DIN |  |  |  |  |
|-------------------------|---|---|--|---|
| ISO Standard | Oxygen | Vacuum | Medical Air | Nitrous Oxide |
| 1/8" NPT |  107916313 3916E-O-NT2 |  107916323 3916E-V-NT2 |  107916333 3916E-A-NT2 |  3916E-2-NT2 |
| 1/4" or 5/16" Hose Barb |  107916210 3916E-O-HB4 |  107916220 3916E-V-HB5 |  107916231 3916E-A-HB4 |  107916240 3916E-2-HB4 |
| Neutral Color | Oxygen | Vacuum | Medical Air | Nitrous Oxide |
| 1/4" or 5/16" Hose Barb |  107916216 3916B-O-HB4 |  107916226 3916B-V-HB5 |  107916231 3916E-A-HB4 | |

Couplers

BS Couplers

| BS |  |  |  |  |
|-------------------------|---|---|--|---|
| ISO Standard | Oxygen | Vacuum | Medical Air | Nitrous Oxide |
| 1/8" NPT |  3917E-O-NT2 |  3917E-V-NT2 |  3917E-A-NT2 |  3917E-2-NT2 |
| 1/4" or 5/16" Hose Barb |  204391710 3917E-O-HB4 |  204391725 3917E-V-HB5 |  204391730 3917E-A-HB4 |  204391740 3917E-2-HB4 |

JIS Couplers

| JIS |  |  |  |  |
|-------------------------|---|---|--|--|
| JIS Standard | Oxygen | Vacuum | Medical Air | Nitrous Oxide |
| 1/8" NPT |  107918012 3918E-O-NT2 |  107918021 3918E-V-NT2 |  107918031 3918E-A-NT2 |  3918E-2-NT2 |
| 1/4" or 5/16" Hose Barb |  107918014 3918E-O-HB4 |  107918020 3918E-V-HB5 |  107918030 3918E-A-HB4 |  3918E-2-HB4 |

Medical Hoses



34-Green



34E-White



34U-Yellow



34E-Blue



34E-Grey



34E-Black

| Part Number | Model Number | Standard | Color | USA Color Service | ISO and GB Color Service |
|-------------|-----------------------------|---------------------------|-------------|-------------------|--------------------------|
| 304460401 | Reel-34-Green ID 1/4 | φ 1/4" (φ 6.3 mm), 100 m | Green | Oxygen | / |
| 304460411 | Reel-34E-White ID 1/4 | φ 1/4" (φ 6.3 mm), 100 m | White | / | Oxygen |
| 304460491 | Reel-34E-Yellow ID 1/4 | φ 1/4" (φ 6.3 mm), 100 m | Yellow | Medical Air | / |
| 304460441 | Reel-34E-Blue ID 1/4 | φ 1/4" (φ 6.3 mm), 100 m | Blue | Nitrous Oxide | Nitrous Oxide |
| 304460461 | Reel-34E-Grey ID 1/4 | φ 1/4" (φ 6.3 mm), 100 m | Grey | Carbon Dioxide | Carbon Dioxide |
| 304460501 | Reel-34E-Black ID 1/4 | φ 1/4" (φ 6.3 mm), 100 m | Black | Nitrogen | Nitrogen |
| 304460431 | Reel-34E-Black+White ID 1/4 | φ 1/4" (φ 6.3 mm), 100 m | Black+White | / | Medical Air |
| 304460421 | Reel-35E-Yellow ID 5/16 | φ 5/16" (φ 8 mm), 100 m | Yellow | / | Vacuun |
| 304460451 | Reel-35E-White ID 5/16 | φ 5/16" (φ 8 mm), 100 m | White | Vacuum | / |

Hose Ferrules

304500240 34E-13 (1/4")
304500250 35E-14.2 (5/16")



Hose Detachable Clips

304501020 34K-D12



Hose Clamp

304501070 34K-11-13



Hose Barb

C55100271 FNT2-HB4
C55100251 FNT2-HB5

Suitable for DISS adapter with check valve, or joint which can not be directly used in the machine to lock the pressure of the machine.



90 Degrees Direction Hose

301100050 1-8FNPT-hose 90
90° Convert, 1/8"-FNPT-HB4



Convert

(Male 1/8" NPT to Hose Barb)

301161150 NT2-HB4
1/8"-FNPT, 1/4"

