

USA Gas Console Outlet

90-Degree DISS (Diameter Index Safety System)

- Accepts Ohmeda®, Chemetron® or Puritan-Bennett® quick connect and DISS gas specific adapters
- Indexed to prevent interchangeability of gas services
- Easy conversion of quick connection or DISS type latch valve assemblies
- 360° swivel inlet for easy installation
- 100% hydrostatically tested
- Complies with NFPA 99 and CGA G-4.1 standards
- UL Listed



Ohmeda® Compatible Shown

British Medical Gas Outlets

BSO Series

GENTEC British Standard Outlets, BSO series, are for use with medical gas systems. The outlets are designed with safety in mind and comply fully with EN 737-1:1998, and HTM 2022. Outlets are manufactured in an ISO 9001:2000 and ISO 13485 registered facility.

- Complies with EN 737-1:1998, HTM 2022, BS5682-1998
- Color-coded for easy gas service identification according to ISO32
- 100% tested and individually packaged
- Gas-specific adapters
- Swivel pipe design for easy installation
- 12mm piping
- CE listed



A Console



B Console

German Medical Gas Outlets

GSO Series

GENTEC German Standard Outlets, GSO series, are for use with medical gas systems. The outlets are designed with safety in mind and comply fully with DIN EN 737-1:1998 and DIN 13620-2. Outlets are manufactured in an ISO 9001:2000 and ISO 13485 registered facility.

- Complies with DIN EN 737-1:1998 and DIN 13620-2
- Color-coded for easy gas service identification according to ISO32
- Stops flow of gas before detachment of probe
- Internal gas-specific safety socket
- 100% tested and individually packaged
- 8mm piping
- CE listed



3816E-O

3826E-O

3826E-O-C

3826E-O-W

French Medical Gas Outlets

F50 Series

GENTEC French Standard Outlets, F50 series, are for use with medical gas systems. The outlets are designed with safety in mind and fully comply with NF EN 737-1:1998. Outlets are manufactured in an ISO 9001 and ISO 13485 certified facility.

- Gas specific adapters
- Color coded for easy gas service identification according to ISO32
- Complies with NF EN 737-1:1998 and NF S 90-116
- UL listed
- 100% tested and individually packaged
- CE listed



F50 (Afnor Style)

Australian Medical Gas Outlets

GENTEC Australian Outlets, SIS series, are used for medical gas pipeline. The outlets are designed with safety in mind and fully comply with EN 737-1:1998. Outlets are manufactured in an ISO 9001 and ISO 13485 certified facility.

- Gas specific adapters
- Color coded for easy gas service identification according to ISO32
- Complies with EN 737-1:1998
- Designed in accordance to S.I.S. (Sleeve Indexed Safety System as per AS2896)
- 100% tested and individually packaged

