


Medical Hose Assemblies

Features

- Designed for centralized medical gas distribution in hospitals, operating rooms, dental offices and clinics
- Maximum working pressure: 200 psi
- Connectors made from brass or stainless steel
- Available in custom lengths and configurations add the length required with “-xx”

Example:

Model number: 34U-OXY-C-OH-DH-12. For a 1/4" green medical hose with ohmeda adapter (male) and DISS Handwheel adapter (male), Oxygen with a 12 inch hose.



Ordering Information

EX: 34	U -	OXY -	C -	OH - DS -	12
Hose	Color Code	Gas Code	Conductive / Non-Conductive	Inlet / Outlet Connection	Hose Length in inches
34: 1/4"	U: US NFPA 99	OXY: Oxygen	C: Conductive	OH: Ohmeda® adapter (male)	6 (0.5 ft)
35: 5/16"	E: ISO standard	VAC: Vacuum	NC: Non-Conductive	DS: DISS hex adapter (male)	12 (1 ft)
		AIR: Medical Air		DH: DISS Handwheel adapter (male)	18 (1.5 ft)
		N2O: Nitrous Oxide		NS: NIST adapter	24 (2 ft)
		WAGD: WAGD		CH: Chemetron® adapter (male)	30 (2.5 ft)
		CO2: Carbon Dioxide		PB: Puritan-Bennett® adapter (male)	36 (3 ft)
		N2: Nitrogen		FS: French (NF S 90-116) adapter (male)	42 (3.5 ft)
		A8: Driving Tool Air		GS: German (DIN 13260-2) adapter (male)	48 (4 ft)
				BS: British (BS5682-1998) adapter (male)	54 (4.5 ft)
				JIS: Japanese Style adapter (male)	60 (5 ft)
				SIS: Australian (AS2896) adapter (male)	66 (5.5 ft)
				FOH: Ohmeda® coupler (female)	72 (6 ft)
				FDS: DISS hex coupler (female)	78 (6.5 ft)
				FDH: DISS Handwheel coupler (female)	84 (7 ft)
				FGS: German (DIN 13260-2) coupler (female)	90 (7.5 ft)
				FBS: British (BS5682-1998) coupler (female)	96 (8 ft)
				FJIS: Japanese Style coupler (female)	228 (19 ft)
				FSIS: Australian (AS2896) coupler (female)	234 (19.5 ft)
				SC: Schrader	240 (20 ft)
				FSC: Schrader coupler	
				FNT2: 1/8" female NPT	
				FNT4: 1/4" female NPT	
				NT2: 1/8" male NPT	
				NT4: 1/4" male NPT	
				OE: medstar	
				OQ: Oxyquip	