

Y-Block Assemblies



CHEMETRON®



OHMEDA®



DISS

Features

- Available in Ohmeda®, Chemetron® and DISS
- Can be used with Vacuum Regulators, Flowmeters, Adapters and DISS connections
- Offered in U.S. and ISO color codes
- Assemblies come standard with 1/8" FNPT
- 5-year warranty

Part Numbers

	Ohmeda®	Chemetron®	DISS
OXYGEN	Y-OXY-CPLR-OH	Y-OXY-CPLR-CH	Y-OXY-CHK-DH
MEDICAL AIR	Y-AIR-CPLR-OH	Y-AIR-CPLR-CH	Y-AIR-CHK-DH
VACUUM	Y-VAC-CPLR-OH	Y-VAC-CPLR-CH	Y-VAC-CHK-DH
N₂O	Y-N ₂ O-CPLR-OH	Y-N ₂ O-CPLR-CH	Y-N ₂ O-CHK-DH
WAGD	Y-WAG-CPLR-OH	Y-WAG-CPLR-CH	Y-WAG-CHK-DH
CO₂	Y-CO ₂ -CPLR-OH	Y-CO ₂ -CPLR-CH	Y-CO ₂ -CHK-DH
NITROGEN	N/A	N/A	Y-NIT-CHK-DH